

alpha 1 Spectrin Rabbit mAb

ABZ3812


Product Overview

Name	alpha 1 Spectrin Rabbit mAb
Accession(Primary)	P02549
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human
Calculated MW	Calculated MW: 280 kDa; Observed MW: 280 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC-P-1:200-1:500


Target information (P02549)

Gene Name SPTA1

About Gene Symbol



Copyright© 2024 Zen-bioscience Limited.



Western blot analysis of SPTA1 in K562 lysates
using alpha 1 Spectrin antibody.