

Desmoplakin Rabbit mAb


ABZ3815

Product Overview

Name	Desmoplakin Rabbit mAb
Accession(Primary)	P15924
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Rat
Calculated MW	Calculated MW: 332 kDa; Observed MW: 270 kDa
Recommended Dilution	WB-1:1000-1:5000; ICC/IF-1:500-1:1000; FC-1:100-1:200


Target information (P15924)

Gene Name	DSP
About Gene Symbol	DSP (Desmoplakin) is a protein kinase and signaling regulator involved in adaptive immunity. It is an important research target in cancer, cell signaling, and therapeutic target research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Copyright© 2024 Zen-bioscience Limited.

Western blot analysis of Desmoplakin in A431
lysates using Desmoplakin antibody.



Western blot analysis of Desmoplakin in A431
lysates using Desmoplakin antibody.