


Adenosine Deaminase Rabbit mAb

ABZ3816

Product Overview


Name	Adenosine Deaminase Rabbit mAb
Accession(Primary)	P00813
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Rat
Calculated MW	Calculated MW: 41 kDa; Observed MW: 41 kDa
Recommended Dilution	WB-1:1000-1:5000

Target information (P00813)

Gene Name ADA**About Gene Symbol**

Western blot analysis of ADA in Jurkat lysates
using ADA antibody.

Copyright© 2024 Zen-bioscience Limited.



Western blot analysis of ADA in Jurkat lysates
using Adenosine Deaminase antibody.