

EIF3E Rabbit mAb

ABZ3817


Product Overview

Name	EIF3E Rabbit mAb
Accession(Primary)	P60228
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 52 kDa; Observed MW: 45 kDa
Recommended Dilution	WB-1:1000-1:5000

Target information (P60228)


Gene Name EIF3E

About Gene Symbol



Western blot analysis of eIF3e in Jurkat, 293T lysates using eIF3e antibody.

Copyright© 2024 Zen-bioscience Limited.



Western blot analysis of eIF3e in (1) 293T lysates;
(2) Jurkat lysates using eIF3E antibody.