

Slit2 Rabbit mAb

ABZ3818


Product Overview

Name	Slit2 Rabbit mAb
Accession(Primary)	O94813
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 170 kDa;Observed MW:58kDa
Recommended Dilution	WB-1:1000-1:5000; IHC-P-1:100-1:200

Target information (O94813)


Gene Name SLIT2

About Gene Symbol



Western blot analysis of Slit2 in SH-SY5Y lysates using Slit2 antibody.

Copyright© 2024 Zen-bioscience Limited.



Western blot analysis of Slit2 in SH-SY5Y lysates using Slit2 antibody.