

COUP TF1 Rabbit mAb


ABZ3821

Product Overview

Name	COUP TF1 Rabbit mAb
Accession(Primary)	P10589
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 46 kDa; Observed MW: 46 kDa
Recommended Dilution	WB-1:500-1:1000; IHC-P-1:100-1:200; ICC/IF-1:100-1:200; FC-1:200-1:500


Target information (P10589)

Gene Name	NR2F1
About Gene Symbol	




Copyright© 2024 Zen-bioscience Limited.

Western blot analysis of COUP TF1 in 293 lysates using COUP TF1 antibody.



Western blot analysis of COUP TF1 in HEK293 lysates using COUP TF1 antibody.



Western blot analysis of COUP TF1 in 293 lysates using COUP TF1 antibody.