

ERG Rabbit mAb


ABZ3822

Product Overview

Name	ERG Rabbit mAb
Accession(Primary)	P11308
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 54 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC-P-1:200-1:1000; ICC/IF-1:200-1:1000; FC-1:1000-1:5000


Target information (P11308)

Gene Name	ERG
About Gene Symbol	




Copyright© 2024 Zen-bioscience Limited.

Western blot analysis of ERG in Jurkat lysates
using ERG antibody.



Western blot analysis of ERG in HeLa lysates
using ERG antibody.



Immunohistochemistry analysis of paraffin-embedded Human kidney using ERG antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.