

NGFR Rabbit mAb


ABZ3823

Product Overview

Name	NGFR Rabbit mAb
Accession(Primary)	P08138
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 45 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC-P-1:100-1:200; ICC/IF-1:20-1:50; FC-1:50-1:100; IP-1:20-1:50


Target information (P08138)

Gene Name	NGFR
About Gene Symbol	




Copyright© 2024 Zen-bioscience Limited.

Western blot analysis of p75 NGF Receptor in PC12 lysates using p75 NGF Receptor antibody.



Western blot analysis of NGFR in (1) C6 lysates; (2) PC-12 lysates using NGFR antibody.



Immunohistochemistry analysis of paraffin-embedded Human glioma using NGFR antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.