

UQCRC2 Rabbit mAb


ABZ3825

Product Overview

Name	UQCRC2 Rabbit mAb
Accession(Primary)	P22695
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 48 kDa; Observed MW: 48 kDa
Recommended Dilution	WB-1:500-1:1000; IHC-P-1:100-1:200; ICC/IF-1:100-1:200; FC-1:50-1:100; IP-1:20-1:50


Target information (P22695)

Gene Name	UQCRC2
About Gene Symbol	UQCRC2 (Ubiquinol-Cytochrome C Reductase Core Protein 2) is a metabolic enzyme involved in cellular metabolism. It is an important research target in metabolic disease, cancer metabolism, and drug discovery and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Copyright© 2024 Zen-bioscience Limited.

Western blot analysis of Complex III Subunit 2 in HepG2, PC-12 lysates using Complex III Subunit 2 antibody.



Western blot analysis of Complex III Subunit 2 in HEK293, Rat brain lysates using Complex III Subunit 2 antibody.