

## MMP2 Rabbit mAb


ABZ3826

### Product Overview


|                             |                                                                         |
|-----------------------------|-------------------------------------------------------------------------|
| <b>Name</b>                 | MMP2 Rabbit mAb                                                         |
| <b>Accession(Primary)</b>   | P08253                                                                  |
| <b>Host Species</b>         | Rabbit                                                                  |
| <b>Clonality</b>            | Monoclonal                                                              |
| <b>Isotype</b>              | IgG                                                                     |
| <b>Species Reactivity</b>   | Human, Mouse, Rat                                                       |
| <b>Calculated MW</b>        | Calculated MW: 74 kDa; Observed MW: 74 kDa                              |
| <b>Recommended Dilution</b> | WB-1:1000-1:5000; ICC/IF-1:100-1:200;<br>FC-1:200-1:500; IHC-1:50-1:100 |

### Target information (P08253)

|                          |      |
|--------------------------|------|
| <b>Gene Name</b>         | MMP2 |
| <b>About Gene Symbol</b> |      |



Western blot analysis of MMP2 in (1)L6 lysates;(2)HeLa lysates using MMP2 antibody.



Immunohistochemistry analysis of paraffin-embedded Human breast adenocarcinoma using MMP2 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.