

HMGA1 Rabbit mAb


ABZ3827

Product Overview

Name	HMGA1 Rabbit mAb
Accession(Primary)	P17096
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 12 kDa; Observed MW: 17 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC-P-1:200-1:500; ICC/IF-1:200-1:500; FC-1:50-1:100

Target information (P17096)

Gene Name	HMGA1
About Gene Symbol	



Western blot analysis of HMGA1 in SK-OV-3
lysates using HMGA1 antibody.

Copyright© 2024 Zen-bioscience Limited.