

Myogenin Rabbit mAb


ABZ3828

Product Overview

Name	Myogenin Rabbit mAb
Accession(Primary)	P15173
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 25 kDa; Observed MW: 36 kDa
Recommended Dilution	WB-1:500-1:1000; IHC-P-1:200-1:500; ICC/IF-1:100-1:1000

Target information (P15173)

Gene Name	MYOG
About Gene Symbol	



Western blot analysis of Myogenin in C2C12
lysates using Myogenin antibody.