


N Cadherin Rabbit mAb

ABZ3831


Product Overview

Name	N Cadherin Rabbit mAb
Accession(Primary)	P19022
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 100 kDa; Observed MW: 140 kDa
Recommended Dilution	WB-1:2000-1:20000; IHC-1:200-1:1000


Target information (P19022)

Gene Name CDH2**About Gene Symbol**

Copyright© 2024 Zen-bioscience Limited.



Western blot analysis of NCadherin in HeLa lysates using N Cadherin antibody.



Immunohistochemistry analysis of paraffin-embedded Human colon using N Cadherin antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.