

Brd4 Rabbit mAb


ABZ3832

Product Overview


Name	Brd4 Rabbit mAb
Accession(Primary)	O60885
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 152 kDa; Observed MW: 152 kDa
Recommended Dilution	WB-1:2000-1:20000; IHC-1:50-1:200; ICC/IF-1:100-1:200; FC-1:20-1:50; IP-1:20

Target information (O60885)


Gene Name	BRD4
About Gene Symbol	




Western blot analysis of BRD4 in HeLa lysates
using BRD4 antibody.



Western blot analysis of Brd4 in HeLa lysates
using BRD4 antibody.



Immunohistochemistry analysis of paraffin-embedded Human kidney using Brd4 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.



Immunofluorescence analysis of BRD4 in HeLa using Brd4 antibody.