

AMPK beta 1 (1A7) Mouse mAb


ABZ3833

Product Overview


Name	AMPK beta 1 (1A7) Mouse mAb
Accession(Primary)	Q9Y478
Host Species	Mouse
Clonality	Monoclonal
Isotype	IgG2a
Species Reactivity	Human, Mouse, Rat, Cat
Calculated MW	Calculated MW: 30 kDa; Observed MW: 38 kDa
Recommended Dilution	WB-1:500-1:1000; IHC-1:50-1:100; IF-1:50-1:200; IP-1:20

Target information (Q9Y478)


Gene Name	PRKAB1
About Gene Symbol	




Western blot analysis of AMPK beta 1 in 3T3, HeLa, PC-12, COS7 and MDA-MB-468 lysates using AMPK beta 1 antibody.



Immunocytochemistry analysis of AMPK beta 1 (1A7) in HeLa using AMPK beta 1 antibody.



Immunoprecipitation analysis of AMPK beta 1 (1A7) in HeLa lysates using AMPK beta 1 antibody.



Immunohistochemistry analysis of paraffin-embedded Breast cancer using AMPK beta 1 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.