


AMPK beta 1 Mouse mAb

ABZ3834

Product Overview

Name	AMPK beta 1 Mouse mAb
Accession(Primary)	Q9Y478
Host Species	Mouse
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human
Calculated MW	Calculated MW: 30 kDa; Observed MW: 38 kDa
Recommended Dilution	WB-1:500-1:1000

Target information (Q9Y478)

Gene Name PRKAB1**About Gene Symbol**

Western blot analysis of AMPK β 1 in HeLa, COS7 lysates using AMPK beta 1 antibody.

