


TRIM33 (3H7) Mouse mAb

ABZ3835

Product Overview

Name	TRIM33 (3H7) Mouse mAb
Accession(Primary)	Q9UPN9
Host Species	Mouse
Clonality	Monoclonal
Isotype	IgG2a
Species Reactivity	Human
Calculated MW	Calculated MW: 123 kDa; Observed MW: 140 kDa
Recommended Dilution	WB-1:500-1:1000

Target information (Q9UPN9)

Gene Name TRIM33**About Gene Symbol**

Western blot analysis of TIF1 gamma / TRIM33 in MCF-7 and HeLa lysates using TIF1 gamma / TRIM33 antibody.

