

Hydroxyacid Oxidase 1 (4B4) Mouse mAb

ABZ3836


Product Overview

Name	Hydroxyacid Oxidase 1 (4B4) Mouse mAb
Accession(Primary)	Q9UJM8
Host Species	Mouse
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Mouse, Rat
Calculated MW	Calculated MW: 41 kDa; Observed MW: 41 kDa
Recommended Dilution	WB-1:500-1:1000; IHC-1:50-1:100; IF-1:50-1:200


Target information (Q9UJM8)

Gene Name HAO1


About Gene Symbol




Immunohistochemistry analysis of paraffin-embedded Human liver tissue using Hydroxyacid Oxidase 1 (4B4) antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval. Negative control was used by secondary antibody only.




Immunohistochemistry analysis of paraffin-embedded Human tonsils using Hydroxyacid Oxidase 1 (4B4) antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval. Negative control was used by secondary antibody only.




Immunohistochemistry analysis of paraffin-embedded mouse heart tissue using HAO1 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval. Negative control was used by secondary antibody only.



Immunofluorescence analysis of Hydroxyacid Oxidase 1 (4B4) in Human appendix using Hydroxyacid Oxidase 1 (4B4) antibody (Mix) (red), and DAPI (blue).



Immunofluorescence analysis of Hydroxyacid Oxidase 1 in mouse spleen using HAO1 antibody (Mix) (red), and DAPI (blue).



Western blot analysis of Hydroxyacid Oxidase 1 in mouse Liver, rat Liver lysates using HAO1 antibody.