


AMACR (2F3) Mouse mAb

ABZ3837


Product Overview

Name	AMACR (2F3) Mouse mAb
Accession(Primary)	Q9UHK6
Host Species	Mouse
Clonality	Monoclonal
Isotype	IgG1
Species Reactivity	Human, Rat
Calculated MW	Calculated MW: 42 kDa; Observed MW: 42 kDa
Recommended Dilution	WB-1:500-1:1000; IF-1:50-1:200

Target information (Q9UHK6)

Gene Name AMACR**About Gene Symbol**

Western blot analysis of AMACR(Cterminus) in rat Liver lysates using AMACR(Cterminus) antibody.



Immunocytochemistry analysis of AMACR (2F3) in Jurkat using AMACR (Cterminus) antibody.