

## SWAP70 (2E9) Mouse mAb


ABZ3838

### Product Overview


|                             |                                            |
|-----------------------------|--------------------------------------------|
| <b>Name</b>                 | SWAP70 (2E9) Mouse mAb                     |
| <b>Accession(Primary)</b>   | Q9UH65                                     |
| <b>Host Species</b>         | Mouse                                      |
| <b>Clonality</b>            | Monoclonal                                 |
| <b>Isotype</b>              | IgG2b                                      |
| <b>Species Reactivity</b>   | Human, Mouse, Rat, Cat                     |
| <b>Calculated MW</b>        | Calculated MW: 69 kDa; Observed MW: 69 kDa |
| <b>Recommended Dilution</b> | WB-1:500-1:1000; IP-1:20                   |

### Target information (Q9UH65)

|                          |                                                                                                                                                                                                                                                                                                                    |
|--------------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Gene Name</b>         | SWAP70                                                                                                                                                                                                                                                                                                             |
| <b>About Gene Symbol</b> | SWAP70 (Switching B Cell Complex Subunit SWAP70) is a protein kinase and signaling regulator involved in adaptive immunity. It is an important research target in cancer, cell signaling, and therapeutic target research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets. |



Western blot analysis of SWAP70 in COS7, 293T, C6, 3T3 and HeLa lysates using SWAP70 antibody.



Immunoprecipitation analysis of Switch Associated Protein 70 (2E9) in HeLa lysates using SWAP70 antibody.