

Lin28A (2C1) Mouse mAb

ABZ3840


Product Overview

Name	Lin28A (2C1) Mouse mAb
Accession(Primary)	Q9H9Z2
Host Species	Mouse
Clonality	Monoclonal
Isotype	IgG2a
Species Reactivity	Human, Mouse
Calculated MW	Calculated MW: 23 kDa; Observed MW: 26 kDa
Recommended Dilution	WB-1:500-1:1000


Target information (Q9H9Z2)

Gene Name LIN28A

About Gene Symbol



Western blot analysis of LIN28 in NTERA and F9 lysates using LIN28A antibody.



Western blot analysis of LIN28A (2C1) in LIN28B and LIN28A recombinant antigen fragments using LIN28A antibody.