

Lin28A (9G2) Mouse mAb

ABZ3841


Product Overview

Name	Lin28A (9G2) Mouse mAb
Accession(Primary)	Q9H9Z2
Host Species	Mouse
Clonality	Monoclonal
Isotype	IgG1
Species Reactivity	Human
Calculated MW	Calculated MW: 23 kDa; Observed MW: 23 kDa
Recommended Dilution	WB-1:500-1:1000; IF-1:50-1:200


Target information (Q9H9Z2)

Gene Name LIN28A


About Gene Symbol



Merge




LIN28




DAPI

Immunofluorescence analysis of LIN28A (9G2) in Hela cells were transfected with pMX construct of Human LIN28 using LIN28A (9G2) antibody



Immunofluorescence analysis of LIN28A (9G2) in NTERA-2 cells using LIN28A (9G2) antibody (green). Blue: DRAQ5 fluorescent DNA dye.



Western blot analysis of LIN28A (9G2) in NTERA-2 lysates using LIN28A (9G2) antibody.