


Wnt5b (3D1) Mouse mAb

ABZ3845

Product Overview

Name	Wnt5b (3D1) Mouse mAb
Accession(Primary)	Q9H1J7
Host Species	Mouse
Clonality	Monoclonal
Isotype	IgG2b
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 40 kDa; Observed MW: 45 kDa
Recommended Dilution	WB-1:500-1:1000

Target information (Q9H1J7)

Gene Name WNT5B**About Gene Symbol**

Western blot analysis of Wnt5b in CHO-K1 lysates (B) and CHO-K1 transfected by Wnt5b (A) lysates using Wnt5b antibody.

