

FOXP3 Mouse (7H9) mAb

ABZ3847


Product Overview


Name	FOXP3 Mouse (7H9) mAb
Accession(Primary)	Q9BZS1
Host Species	Mouse
Clonality	Monoclonal
Isotype	IgG1
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 47 kDa; Observed MW: 47 kDa
Recommended Dilution	WB-1:500-1:1000; IF-1:50-1:200

Target information (Q9BZS1)

Gene Name FOXP3

About Gene Symbol





Immunocytochemistry analysis of FOXP3 mouse (7H9) mAbCHO-K1 which transfected by FOXP3 using FOXP3 antibody.