

TBLR1/TBL1XR1 (4F3) Mouse mAb

ABZ3848


Product Overview

Name	TBLR1/TBL1XR1 (4F3) Mouse mAb
Accession(Primary)	Q9BZK7
Host Species	Mouse
Clonality	Monoclonal
Isotype	IgG1
Species Reactivity	Human, Mouse
Calculated MW	Calculated MW: 56 kDa; Observed MW: 60 kDa
Recommended Dilution	WB-1:500-1:1000; IHC-1:50-1:100; IF-1:50-1:200


Target information (Q9BZK7)

Gene Name TBL1XR1


About Gene Symbol




Western blot analysis of TBLR1 in mouse brain, mouse heart and K562 lysates using TBLR1 antibody.



Immunocytochemistry analysis of TBL1XR1 (4F3) in HeLa using TBLR1 antibody.



Immunohistochemistry analysis of paraffin-embedded human breast cancer using TBLR1 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.



Western blot analysis of TBL1XR1 (4F3) in HeLa, K562 lysates using TBL1XR1 (4F3) antibody

