

TBLR1/TBL1XR1 Mouse mAb

ABZ3849

Product Overview

Name	TBLR1/TBL1XR1 Mouse mAb
Accession(Primary)	Q9BZK7
Host Species	Mouse
Clonality	Monoclonal
Isotype	IgG1
Species Reactivity	Human, Cat
Calculated MW	Calculated MW: 56 kDa; Observed MW: 60 kDa
Recommended Dilution	WB-1:500-1:1000; IHC-1:50-1:100; IF-1:50-1:200

Target information (Q9BZK7)

Gene Name TBL1XR1

About Gene Symbol



