

EFHD1 (5C10) Mouse mAb

ABZ3850


Product Overview

Name	EFHD1 (5C10) Mouse mAb
Accession(Primary)	Q9BUP0
Host Species	Mouse
Clonality	Monoclonal
Isotype	IgG1
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 27 kDa; Observed MW: 27 kDa
Recommended Dilution	WB-1:500-1:1000; IHC-1:50-1:100; IF-1:50-1:200


Target information (Q9BUP0)

Gene Name EFHD1


About Gene Symbol




Immunohistochemistry analysis of paraffin-embedded Human uterus tissue using EFHD1 (5C10) antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval. Negative control was used by secondary antibody only.




Immunohistochemistry analysis of paraffin-embedded Human tonsils using EFHD1 (5C10) antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval. Negative control was used by secondary antibody only.




Immunohistochemistry analysis of paraffin-embedded mouse testis tissue using EFHD1 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval. Negative control was used by secondary antibody only.



Immunofluorescence analysis of EFHD1 (5C10) in mouse lung using EFHD1 (5C10) antibody (3G2) (red), and DAPI (blue).



Western blot analysis of EFHD1 in mouse spleen, rat spleen tissue lysates using EFHD1 (5C10) antibody.



Immunofluorescence analysis of EFHD1 in HeLa using EFHD1 (5C10) antibody(Left) and DAPI (Right).