


RAD9A Mouse mAb

ABZ3852


Product Overview

Name	RAD9A Mouse mAb
Accession(Primary)	Q99638
Host Species	Mouse
Clonality	Monoclonal
Isotype	IgG2b
Species Reactivity	Human, Mouse, Rat, Cat
Calculated MW	Calculated MW: 43 kDa; Observed MW: 60 kDa
Recommended Dilution	WB-1:500-1:1000; IP-1:20

Target information (Q99638)

Gene Name RAD9A**About Gene Symbol**

Western blot analysis of RAD9A in HeLa, MCF-7, 3T3, COS7 and rat testis lysates using RAD9A antibody.



Immunoprecipitation analysis of Rad9 in HeLa and 3T3 lysates using RAD9A antibody.