


CDX2 (4B9) Mouse mAb

ABZ3853


Product Overview

Name	CDX2 (4B9) Mouse mAb
Accession(Primary)	Q99626
Host Species	Mouse
Clonality	Monoclonal
Isotype	IgG2b
Species Reactivity	Human, Mouse, Cat
Calculated MW	Calculated MW: 34 kDa; Observed MW: 34 kDa
Recommended Dilution	WB-1:500-1:1000; IF-1:50-1:200; IP-1:20


Target information (Q99626)

Gene Name CDX2**About Gene Symbol**

Western blot analysis of CDX2 in HeLa, Caco2, C2C12, SW480, COS7, HepG2 and 3T3 lysates using CDX2 antibody.



Immunocytochemistry analysis of CDX2 (4B9) in HeLa using CDX2 antibody.



Immunoprecipitation analysis of CDX2 (4B9) in HeLa cell lysates using CDX2 antibody.