

CDX2 (4F12) Mouse mAb

ABZ3854


Product Overview

Name	CDX2 (4F12) Mouse mAb
Accession(Primary)	Q99626
Host Species	Mouse
Clonality	Monoclonal
Isotype	IgG1
Species Reactivity	Human, Cat
Calculated MW	Calculated MW: 34 kDa; Observed MW: 34 kDa
Recommended Dilution	WB-1:500-1:1000; IF-1:50-1:200; IP-1:20


Target information (Q99626)

Gene Name CDX2


About Gene Symbol



Western blot analysis of CDX2 in HeLa, COS7 and K562 lysates using CDX2 antibody.



Immunocytochemistry analysis of CDX2 (4F12) in HeLa using CDX2 antibody.



Immunoprecipitation analysis of CDX2 (4F12) in HeLa lysates using CDX2 antibody.