


UHRF1 (1H9) Mouse mAb

ABZ3855

Product Overview

Name	UHRF1 (1H9) Mouse mAb
Accession(Primary)	Q96T88
Host Species	Mouse
Clonality	Monoclonal
Isotype	IgG2b
Species Reactivity	Human, Mouse
Calculated MW	Calculated MW: 90 kDa; Observed MW: 90 kDa
Recommended Dilution	WB-1:500-1:1000

Target information (Q96T88)

Gene Name UHRF1**About Gene Symbol**

Western blot analysis of UHRF1 in mouse heart, 293T, HeLaNE and HeLa lysates using UHRF1 antibody.

