

PRMT6 (2C3) Mouse mAb


ABZ3856

Product Overview


Name	PRMT6 (2C3) Mouse mAb
Accession(Primary)	Q96LA8
Host Species	Mouse
Clonality	Monoclonal
Isotype	IgG1
Species Reactivity	Human
Calculated MW	Calculated MW: 42 kDa; Observed MW: 42 kDa
Recommended Dilution	WB-1:500-1:1000; IHC-1:50-1:100; IF-1:50-1:200

Target information (Q96LA8)


Gene Name	PRMT6
About Gene Symbol	PRMT6 (Protein Arginine Methyltransferase 6) is a methyltransferase involved in chromatin remodeling and epigenetics. It is an important research target in metabolic disease, cancer metabolism, and drug discovery and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Western blot analysis of PRMT6 in HeLa and 293T lysates using PRMT6 antibody.



Immunofluorescence analysis of PRMT6 (2C3) in HeLa using PRMT6 antibody.



Immunohistochemistry analysis of paraffin-embedded Human breast cancer using PRMT6 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.