

## PRMT6 Mouse mAb


ABZ3857

### Product Overview

|                             |                                                |
|-----------------------------|------------------------------------------------|
| <b>Name</b>                 | PRMT6 Mouse mAb                                |
| <b>Accession(Primary)</b>   | Q96LA8                                         |
| <b>Host Species</b>         | Mouse                                          |
| <b>Clonality</b>            | Monoclonal                                     |
| <b>Isotype</b>              | IgG1                                           |
| <b>Species Reactivity</b>   | Human, Cat                                     |
| <b>Calculated MW</b>        | Calculated MW: 42 kDa; Observed MW: 42 kDa     |
| <b>Recommended Dilution</b> | WB-1:500-1:1000; IHC-1:50-1:100; IF-1:50-1:200 |

### Target information (Q96LA8)

|                          |                                                                                                                                                                                                                                                                                                               |
|--------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Gene Name</b>         | PRMT6                                                                                                                                                                                                                                                                                                         |
| <b>About Gene Symbol</b> | PRMT6 (Protein Arginine Methyltransferase 6) is a methyltransferase involved in chromatin remodeling and epigenetics. It is an important research target in metabolic disease, cancer metabolism, and drug discovery and is widely studied to elucidate cellular mechanisms and identify therapeutic targets. |



Western blot analysis of PRMT6 in COS7, Jurkat, C6, 3T3 and HeLa lysates using PRMT6 antibody.