

SIRT1 (2G10) Mouse mAb

ABZ3858


Product Overview

Name	SIRT1 (2G10) Mouse mAb
Accession(Primary)	Q96EB6
Host Species	Mouse
Clonality	Monoclonal
Isotype	IgG2b
Species Reactivity	Human
Calculated MW	Calculated MW: 82 kDa; Observed MW: 120 kDa
Recommended Dilution	WB-1:500-1:1000; IP-1:20


Target information (Q96EB6)

Gene Name SIRT1

About Gene Symbol



Western blot analysis of SIRT1 in HeLa, Jurkat, Ramos, A549 and CEM lysates using SIRT1 antibody.



Immunoprecipitation analysis of SIRT1 (2G10) in HeLa lysates using SIRT1 antibody.