

USP7 (3E4) Mouse mAb


ABZ3860

Product Overview

Name	USP7 (3E4) Mouse mAb
Accession(Primary)	Q93009
Host Species	Mouse
Clonality	Monoclonal
Isotype	IgG1
Species Reactivity	Human
Calculated MW	Calculated MW: 128 kDa; Observed MW: 135 kDa
Recommended Dilution	WB-1:500-1:1000

Target information (Q93009)

Gene Name	USP7
About Gene Symbol	USP7 (Ubiquitin Specific Peptidase 7) is a deubiquitinase involved in ubiquitin-proteasome pathway. It is an important research target in infectious disease and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Western blot analysis of HAUSP / USP7 in HeLa and K562 lysates using HAUSP / USP7 antibody.