

SMARCC1/BAF155 Mouse mAb


ABZ3861

Product Overview

Name	SMARCC1/BAF155 Mouse mAb
Accession(Primary)	Q92922
Host Species	Mouse
Clonality	Monoclonal
Isotype	IgG1
Species Reactivity	Human
Calculated MW	Calculated MW: 123 kDa; Observed MW: 155 kDa
Recommended Dilution	WB-1:500-1:1000

Target information (Q92922)

Gene Name	SMARCC1
About Gene Symbol	SMARCC1 (SWI/SNF Complex Subunit SMARCC1) is a structural cytoskeletal protein involved in chromatin remodeling and epigenetics. It is an important research target in cell biology, cancer metastasis, and cytoskeletal research and is widely studied to eluc



Western blot analysis of SMARCC1 in K562, Jurkat and HeLa lysates using SMARCC1 antibody.