


GATA6 (3H3) Mouse mAb

ABZ3862

Product Overview

Name	GATA6 (3H3) Mouse mAb
Accession(Primary)	Q92908
Host Species	Mouse
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human
Calculated MW	Calculated MW: 60 kDa; Observed MW: 60 kDa
Recommended Dilution	WB-1:500-1:1000

Target information (Q92908)

Gene Name GATA6**About Gene Symbol**

Western blot analysis of GATA6 in HeLa and 293 lysates using GATA6 antibody.

