

HDAC2 (2D9) Mouse mAb


ABZ3863

Product Overview

Name	HDAC2 (2D9) Mouse mAb
Accession(Primary)	Q92769
Host Species	Mouse
Clonality	Monoclonal
Isotype	IgG2b
Species Reactivity	Human, Mouse, Rat, Cat
Calculated MW	Calculated MW: 55 kDa; Observed MW: 60 kDa
Recommended Dilution	WB-1:500-1:1000; IF-1:50-1:200

Target information (Q92769)

Gene Name	HDAC2
About Gene Symbol	HDAC2 (Histone Deacetylase 2) is a histone deacetylase involved in chromatin remodeling and epigenetics. It is an important research target in epigenetic regulation, chromatin biology, and cancer research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Western blot analysis of HDAC2 in THP-1, SH-SY5Y, K562, COS7, C6, 3T3 and HeLa lysates using HDAC2 antibody.