

MLKL (5E3) Mouse mAb


ABZ3864

Product Overview


Name	MLKL (5E3) Mouse mAb
Accession(Primary)	Q8NB16
Host Species	Mouse
Clonality	Monoclonal
Isotype	IgG1
Species Reactivity	Human
Recommended Dilution	IHC-1:50-1:100

Target information (Q8NB16)

Gene Name	MLKL
About Gene Symbol	MLKL (Mixed Lineage Kinase Domain-Like Protein) is a protein kinase and signaling regulator involved in cell signaling and proliferation. It is an important research target in cancer, cell signaling, and therapeutic target research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Immunohistochemistry analysis of paraffin-embedded Human Tonsil Tissue using MLKL (5E3) antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.



Immunohistochemical analysis of paraffin-embedded Human tonsils using MLKL (5E3) antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.