

SMYD1 (2E11) Mouse mAb

ABZ3865


Product Overview

Name	SMYD1 (2E11) Mouse mAb
Accession(Primary)	Q8NB12
Host Species	Mouse
Clonality	Monoclonal
Isotype	IgG2b
Species Reactivity	Rat
Calculated MW	Calculated MW: 57 kDa; Observed MW: 57 kDa
Recommended Dilution	WB-1:500-1:1000

Target information (Q8NB12)

Gene Name SMYD1

About Gene Symbol



Western blot analysis of SMYD1 in rat skeletal muscle lysates using SMYD1 antibody.

