

SIRT6 (6C9) Mouse mAb

ABZ3866


Product Overview

Name	SIRT6 (6C9) Mouse mAb
Accession(Primary)	Q8N6T7
Host Species	Mouse
Clonality	Monoclonal
Isotype	IgG1
Species Reactivity	Human, Mouse, Rat, Cat
Calculated MW	Calculated MW: 39 kDa; Observed MW: 42 kDa
Recommended Dilution	WB-1:500-1:1000; IF-1:50-1:200; IP-1:20


Target information (Q8N6T7)

Gene Name SIRT6


About Gene Symbol



Western blot analysis of SIRT6 (6C9) in Jurkat, MCF-7, HCT116, K562, COS7, C2C12, 3T3 and HeLa lysates using SIRT6 antibody.



Immunofluorescence analysis of SIRT6 (6C9) in HeLa using SIRT6 antibody.



Immunoprecipitation analysis of SIRT6 (6C9) in HeLa lysates using SIRT6 antibody.