

DBC 1 (3G4) Mouse mAb

ABZ3867


Product Overview

Name	DBC 1 (3G4) Mouse mAb
Accession(Primary)	Q8N163
Host Species	Mouse
Clonality	Monoclonal
Isotype	IgG1
Species Reactivity	Human, Mouse, Rat, Cat
Calculated MW	Calculated MW: 103 kDa; Observed MW: 130 kDa
Recommended Dilution	WB-1:500-1:1000; IF-1:50-1:200; IP-1:20


Target information (Q8N163)

Gene Name CCAR2


About Gene Symbol



Western blot analysis of DBC1 in HeLa, A549, Jurkat, C6, 3T3 and COS7 lysates using DBC1 antibody.



Immunocytochemistry analysis of DBC 1 in HeLa cells using DBC1 antibody.



Immunoprecipitation analysis of DBC 1 (3G4) in HeLa lysates using DBC1 antibody.