

SIRT2 (1D4) Mouse mAb

ABZ3868


Product Overview

Name	SIRT2 (1D4) Mouse mAb
Accession(Primary)	Q8IXJ6
Host Species	Mouse
Clonality	Monoclonal
Isotype	IgG2b
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 43 kDa; Observed MW: 43 kDa
Recommended Dilution	WB-1:500-1:1000; IF-1:50-1:200


Target information (Q8IXJ6)

Gene Name SIRT2


About Gene Symbol



Western blot analysis of SirT2 in PC-12, rat Brain and mouse Brain lysates using SirT2 antibody.



Immunofluorescence analysis of SIRT2 (1D4) in Hela using SirT2 antibody(green) and Histone H3.1 (Phospho-Ser10) (310045) antibody(red).



Western blot analysis of SIRT2 (1D4) in rat Cerebral Cortex and rat Brain lysates using SirT2 antibody.