

CARM1 (2B9) Mouse mAb

ABZ3869


Product Overview

Name	CARM1 (2B9) Mouse mAb
Accession(Primary)	Q86X55
Host Species	Mouse
Clonality	Monoclonal
Isotype	IgG1
Species Reactivity	Human, Mouse
Calculated MW	Calculated MW: 66 kDa; Observed MW: 63 kDa
Recommended Dilution	WB-1:500-1:1000; IP-1:20


Target information (Q86X55)

Gene Name CARM1

About Gene Symbol



Western blot analysis of PRMT4/CARM1 in HeLa, A431 and K562 lysates using PRMT4/CARM1 antibody.



Immunoprecipitation analysis of CARM1 (2B9) in HeLa lysates using PRMT4/CARM1 antibody.