

LYRIC/AEG1 (2E4) Mouse mAb

ABZ3870


Product Overview


Name	LYRIC/AEG1 (2E4) Mouse mAb
Accession(Primary)	Q86UE4
Host Species	Mouse
Clonality	Monoclonal
Isotype	IgG1
Species Reactivity	Human
Calculated MW	Calculated MW: 64 kDa; Observed MW: 64 kDa
Recommended Dilution	WB-1:500-1:1000; IHC-1:50-1:100; FC-1:50-1:100

Target information (Q86UE4)


Gene Name MTDH

About Gene Symbol





Flow cytometry analysis of Hela cells stained with Metadherin antibody(green) and negative control(purple).



Immunohistochemistry analysis of paraffin-embedded Human Liver tissues using Metadherin antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.