

TORC1/CRTC1 Mouse mAb


ABZ3872

Product Overview


Name	TORC1/CRTC1 Mouse mAb
Accession(Primary)	Q6UUV9
Host Species	Mouse
Clonality	Monoclonal
Isotype	IgG2b
Species Reactivity	Human
Calculated MW	Calculated MW: 67 kDa; Observed MW: 78 kDa
Recommended Dilution	WB-1:500-1:1000; IF-1:50-1:200; IP-1:20; FC-1:50-1:100

Target information (Q6UUV9)


Gene Name	CRTC1
About Gene Symbol	




Western blot analysis of MECT1 / Torc1 in HeLa and Jurkat lysates using MECT1 / Torc1 antibody.



Immunocytochemistry analysis of CRTC1 in HeLa using MECT1 / Torc1 antibody.



Immunoprecipitation analysis of CRTC1 in HeLa cell lysates using MECT1 / Torc1 antibody.



Flow Cytometry analysis of CRTC1 in K562 cells using CRTC1 antibody (red). Blue line histogram represents the isotype control.