

Phostensin (4B12) Mouse mAb


ABZ3873

Product Overview

Name	Phostensin (4B12) Mouse mAb
Accession(Primary)	Q6NYC8
Host Species	Mouse
Clonality	Monoclonal
Isotype	IgG1
Species Reactivity	Human
Calculated MW	Calculated MW: 68 kDa; Observed MW: 110 $\times$ 26 kDa
Recommended Dilution	WB-1:500-1:1000

Target information (Q6NYC8)

Gene Name	PPP1R18
About Gene Symbol	



Western blot analysis of Phostensin in HeLa, Raji, Ramos and Leucocyte lysates using Phostensin antibody.