


JMJD6 (3G5) Mouse mAb

ABZ3874


Product Overview

Name	JMJD6 (3G5) Mouse mAb
Accession(Primary)	Q6NYC1
Host Species	Mouse
Clonality	Monoclonal
Isotype	IgG1
Species Reactivity	Human
Calculated MW	Calculated MW: 46 kDa; Observed MW: 62 kDa
Recommended Dilution	WB-1:500-1:1000; IF-1:50-1:200

Target information (Q6NYC1)

Gene Name JMJD6**About Gene Symbol**

Western blot analysis of JMJD6(Nterminus) in HeLa and HL-60 lysates using JMJD6(Nterminus) antibody.



Immunocytochemistry analysis of JMJD6 in HeLa cells using JMJD6(Nterminus) antibody.