


ABCB5 (8D2) Mouse mAb

ABZ3875


Product Overview

Name	ABCB5 (8D2) Mouse mAb
Accession(Primary)	Q2M3G0
Host Species	Mouse
Clonality	Monoclonal
Isotype	IgG1
Species Reactivity	Human
Calculated MW	Calculated MW: 139 kDa; Observed MW: 90 kDa
Recommended Dilution	WB-1:500-1:1000; IF-1:50-1:200


Target information (Q2M3G0)

Gene Name ABCB5**About Gene Symbol**

Western blot analysis of ABCB5 in K562 lysates using ABCB5 antibody.



Western blot analysis of ABCB5 (8D2) in HepG2, K562 lysates using ABCB5 (8D2) antibody.



Immunocytochemistry analysis of ABCB5 (8D2) (green) in Hela using ABCB5 (8D2) antibody, and DAPI(blue)