

KIF7 (2A7) Mouse mAb

ABZ3876


Product Overview


Name	KIF7 (2A7) Mouse mAb
Accession(Primary)	Q2M1P5
Host Species	Mouse
Clonality	Monoclonal
Isotype	IgG1
Species Reactivity	Human, Mouse, Rat
Recommended Dilution	IHC-1:50-1:100; IF-1:50-1:200

Target information (Q2M1P5)

Gene Name KIF7

About Gene Symbol





Immunohistochemistry analysis of paraffin-embedded mouse Kidney tissue using KIF7 (2A7) antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.